

Dental Treatments

The mild abrasive nature of Navajo Brand® pumice makes it an ideal polishing aid to finish acrylic dentures. It is also used in oral hygiene pastes, prophy pastes, and teeth polishing compounds.

Navajo Brand® pumice is used in a water suspension to obtain a smooth surface during polishing.

Navajo Brand® pumice has several products for this application

Recommended Products

NB 012	12-Micron Average Particle Size
FFFF	90% Finer than 44 microns (325 Mesh)
FFF	80% Finer than 44 microns (325 Mesh)
F	70% Finer than 44 microns (325 Mesh)
#4/0	45% Finer than 44 microns (325 Mesh)
#1/2	0.15 to 0.3 mm (50 to 120 Mesh)
#3	0.2 to 0.6 mm (30 to 80 Mesh)

Additional Information at www.CRMinerals.com

To place an order or obtain additional information, please contact CR Minerals at 505-428-2940, or contact your local distributor.

CR Minerals is a worldwide supplier of pumice products to many diverse markets. It operates a state of the art processing facility in Ohkay Owingeh, New Mexico.

Although the information and suggestions in this publication are believed to be correct, CR Minerals makes no representations or warranties as to the accuracy or completeness of this information.

